

## LLLT Series

# M-Shaped Hairline laser therapy device

650nm

The red light penetrates surface cells to reach the hair follicle roots, stimulating the rich melanin there to specifically absorb this soft laser wavelength. This energy absorption increases blood supply to the entire follicle area, accelerates ATP synthesis, and activates follicle regrowth.

Designed for M-shaped baldness and hairline issues, the headband is comfortable with a skin-friendly lining and easy to use anytime, anywhere. The breathable fabric keeps it dry and snug against your skin, ensuring effective light penetration.



## Feature

### Accelerate circulation, activate hair follicles

Stimulate roots of hair follicles with soft laser energy

- Accelerate blood circulation
- Accelerate ATP formation
- Activate redevelopment of hair follicles.

## I Specification

Wavelength	650nm
Laser specifications	8pcs
Output power	Class 2M, <1mW
Battery capacity	Lithium polymer battery 400mAh
Product size	Φ40/H13 ±1mm
Weight	17 ±1g
Operating mode	Continuous wave (CW)/ Pulse wave (PW)
Operating environment	Temperature 10~40 degrees; Humidity 30~75%
Transportation and storage environment	Temperature -20~40 degrees; Humidity 10~90%

## 【 Cautions 】

- Avoid light exposure to your eyes directly for protect your eyes.
- Transformer is prohibited for use which may cause dangerous event.
- Do not use fast charging for this product, as it may cause danger.
- Please keep clean on the area of laser output, prevent the dirties affect energy transmission when use.
- Please keep it away from the touch by children or pets in order to avoid eating or touch components accidentally.
- Do not use on pregnant women, who with skin disease, malignant tumor, wear pacemaker, initial wound within 72 hours, physical abnormalities and abnormal sensory nerves. Avoid to use it during the skin is sensitive, fragile and in rash or wounds condition. Please stop using it immediately and consult doctor if you feel any discomfort.
- Don' t use it with electronic medical equipment in order to avoid the malfunction or dangerous accidents.
- Please store it in dry, cool place and away from water. Please pay attention on moisture and falling during use.

### I LECC TECHNOLOGY CO., LTD.

2F No.2, Ziqiang 3rd Rd., Zhongli Dist., Taoyuan City, 320023, Taiwan.  
 TEL | 03 4511967 VAT Number | 16708380  
 FAX | 03 4511923 Email | sales@lecc.com.tw

### I SHANGHAI LECC OPTO CO., LTD.

No.1939, Da Ye Rd., Feng-Xian Dist., Shanghai City, 201402, China.  
 TEL | 021 57407400 TAX Number | 913100007630197034  
 FAX | 021 57407403 Email | sales@lecc.com.tw