

LLLT Series

M-Shaped Hairline laser therapy device

650nm

The red light penetrates the surface cells, reaches and stimulates the root of the hair follicle with rich melanin which absorb laser energy of this wavelength specifically. Increasing the blood supply to whole hair follicle area after absorb energy. Accelerating ATP synthesis and activate the regrowth for hair follicles. Be designed as a hat for head full coverage to achieve full range of irradiation.

Designed for M-shaped baldness and hairline issues. Using headband is comfortable and convenient, it's easy to use anytime and anywhere. The fabrics are breathable and dry which fit your skin comfortably, allowing light penetrates your skin effectively.



Feature

- Free to change continuous on and pulse mode.
- Relax the hair follicles on single area gently.
- Only 17g, it's a slim appearance which is easy to carry.
- The indicator of low battery capacity which remind you to charge battery via Type-C.

Specification

Wavelength	650 ±20nm
Product size	Φ40/H13 ±1mm
Product Weight	17 ±1g
Battery Capacity	Lithium polymer battery 400mAh
Laser Safety Classification	Class 1
Maximum Optical Power	0.311mW
Operating mode	continuous wave (CW)/ pulse wave (PW)
Laser standard	IEC 60825-1:2014
Operating environment	Temperature 10~40 degrees; Humidity 30~75%
Transportation and storage environment	Temperature -20~40 degrees; Humidity 10~90%

【 Cautions 】

- Avoid light exposure to your eyes directly for protect your eyes.
- Transformer is prohibited for use which may cause dangerous event.
- Do not use fast charging for this product, as it may cause danger.
- Please keep clean on the area of laser output and scalp, prevent the dirties affect energy transmission when use.
- Please keep it away from the touch by children or pets in order to avoid eating or touch components accidentally.
- Do not use on pregnant women, who with skin disease, malignant tumor, wear pacemaker, initial wound within 72 hours, physical abnormalities and abnormal sensory nerves. Avoid to use it during the skin is sensitive, fragile and in rash or wounds condition. Please stop using it immediately and consult doctor if you feel any discomfort.
- Don' t use it with electronic medical equipment in order to avoid the malfunction or dangerous accidents.
- Please store it in dry, cool place and away from water. Please pay attention on moisture and falling during use.

LECC TECHNOLOGY CO., LTD.

2F No.2, Ziqiang 3rd Rd., Zhongli Dist., Taoyuan City, 320023, Taiwan.
 TEL | 03 4511967 VAT Number | 16708380
 FAX | 03 4511923 Email | sales@lecc.com.tw

SHANGHAI LECC OPTO CO., LTD.

No.1939, Da Ye Rd., Feng-Xian Dist., Shanghai City, 201402, China.
 TEL | 021 57407400 TAX Number | 913100007630197034
 FAX | 021 57407403 Email | sales@lecc.com.tw