

## LLLT Series

# Nasal laser therapy device

650nm

The nasal care device activates nasal collaterals to improve energy in the nasal area, helping to alleviate dark circles and achieve other beneficial effects. Through gentle irradiation of the nasal cavity with a wavelength of 650nm and the infusion of a significant number of photons, the device utilizes the stimulating principles of light on cells. This process penetrates deep into the skin and blood vessels, enabling cells to activate self-repair functions. It relaxes nasal muscles and alleviates fatigue, effectively saying goodbye to dark circles.



## Feature

### Accelerate circulation, relieve Dark Circles

Health care for nose of user, to achieve the following effects

- Smell and aroma conditioning
- Relax nose muscles
- Relieve fatigue
- Accelerate blood circulation
- Relieves Dark Circles

## I Specification

Wavelength	650nm
Laser specifications	2pcs
Output power	Class 2M, <1mW
Working voltage	3-5 V, 40 mA(MAX.)
Product size	H1200 ±5mm
Weight	47 ±5g
Operating environment	Temperature 10~40 degrees; Humidity 30~75%
Transportation and storage environment	Temperature 0~60 degrees; Humidity 10~90%

## 【 Cautions 】

- Avoid light exposure to your eyes directly for protect your eyes.
- Transformer is prohibited for use which may cause dangerous event.
- Do not use fast charging for this product, as it may cause danger.
- Please keep clean on the area of laser output, prevent the dirties affect energy transmission when use.
- Please keep it away from the touch by children or pets in order to avoid eating or touch components accidentally.
- Do not use on pregnant women, who with skin disease, malignant tumor, wear pacemaker, initial wound within 72 hours, physical abnormalities and abnormal sensory nerves. Avoid to use it during the skin is sensitive, fragile and in rash or wounds condition. Please stop using it immediately and consult doctor if you feel any discomfort.
- Don' t use it with electronic medical equipment in order to avoid the malfunction or dangerous accidents.
- Please store it in dry, cool place and away from water. Please pay attention on moisture and falling during use.

### I LECC TECHNOLOGY CO., LTD.

2F No.2, Ziqiang 3rd Rd., Zhongli Dist., Taoyuan City, 320023, Taiwan.  
 TEL | 03 4511967 VAT Number | 16708380  
 FAX | 03 4511923 Email | sales@lecc.com.tw

### I SHANGHAI LECC OPTO CO., LTD.

No.1939, Da Ye Rd., Feng-Xian Dist., Shanghai City, 201402, China.  
 TEL | 021 57407400 TAX Number | 913100007630197034  
 FAX | 021 57407403 Email | sales@lecc.com.tw